




APZEM PROFESSIONAL BACKDRAFT GROSSING STATION (RB Series)

APZEM Professional series Backdraft Grossing Station were designed to provide Sturdy, customizable and affordable workstation for medium to large Labs. The professional series workstation comes with higher suction capacity and additional features that assists in grossing as well as protects you from inhaling harmful Formalin and other harmful fumes while grossing.



| MODEL NO | SIZE (LENGTH) |
|----------|---------------|
| RB400 | 4 FEET |
| RB500 | 5 FEET |
| RB600 | 6 FEET |
| RB700 | 7 FEET |



STANDARD DESIGN FEATURES

- ▶ All Stainless steel Construction made of thicker 16 Guage SS 304.
- ▶ Inbuilt ventilation system with dual wall construction capable of removing all formalin fumes over work surface and from Technician's breathing zone.
- ▶ Deep Sink with Cold water Faucet.
- ▶ Back Draft fume Removal system.
- ▶ Removable perforated worktop.
- ▶ Powerful LED Lighting.
- ▶ C-Fold paper towel Holder.
- ▶ Polyethylene Dissecting Board for Dissecting and for Photography.
- ▶ Magnetic Tool holder.
- ▶ Cleanable Pre Filter to remove fine particles.
- ▶ Worktable Height : 800mm.
- ▶ Additional shelf for storing.
- ▶ Hand rinsing assembly.
- ▶ Perimeter Rinse assembly.
- ▶ High capacity 0.55 HP waste disposal unit.
- ▶ Transparent Lexan air deflector panels.
- ▶ 40mm height adjustable leg for accurate leveling of work surface.
- ▶ Formalin dispensing system.
- ▶ Provided with additional switch and socket.

OPTIONAL FEATURES AND ACCESSORIES

- ▶ Magnifying Lamp fixture for mounting 3 axis adjustable magnifying Lamp.
- ▶ CPU Mounting bracket.
- ▶ Digital camera holder.
- ▶ Safety splash shield.
- ▶ Adjustable Shelving.
- ▶ Pull out writing pad.
- ▶ Organizer Bin.
- ▶ Optional Carbon filter / potassium permanganate filter.
- ▶ LED Monitor and Keyboard mounting bracket.
- ▶ Higher grade, SS 316 Construction for better chemical resistance and greater life.

** Chose any of 48 accessories listed at end..